

CT Pneumonia Triage


is a machine learning based triage software designed for identifying pneumonia findings from non-contrast chest CT images. It can be used by radiologists for earlier notification and worklist prioritization purposes.

CT Pneumonia Triage uses a deep learning-based algorithm for the detection of pneumonia findings from non-contrast chest CT images. The results, available to the PACS worklist in less than one minute, can be used to notify radiologists that the image includes findings suggestive of pneumonia so that a radiologist can review sooner.

Besides general population pneumonia cases, the AI algorithm has also been trained and tested with COVID-19 pneumonia cases, its high sensitivity (above 90%) indicates it may be used to alert radiologists of suspected pneumonia findings that could be associated with COVID-19. While COVID-19 remains to be a pandemic, this tool may allow for earlier notification of COVID-19-related findings in otherwise asymptomatic patients.

Clinical Workflow

CT Pneumonia Triage can be seamlessly integrated with worklist applications and PACS/VNA (Vendor Neutral Archives).



What it does:

The software may be activated via an application programming interface (API) from a worklist application. Using deep learning techniques trained from thousands of annotated CT images, the software analyzes chest CT images and sends the binary prediction results (suspected pneumonia or not) to the worklist. The suspected pneumonia cases may be highlighted by the worklist application to alert radiologists for earlier treatment.

Patient Name	MRN	Modality	Procedure	Arrival Time
Glazer, Aaron	87634538	X-Ray	CXR ABD	3 hr ago
Miller, Sophia	28394751	X-Ray	CXR Chest	30 min ago
Vloet, Sanne	98726534	CT	CT Chest W/O Contrast	20 min ago
Rice, Jessica	12361524	CT	CT Head W/O Contrast	15 min ago
Li, Angel	26378940	CT	CT Chest W/O Contrast	9 min ago

Standard worklist



Patient Name	MRN	Modality	Procedure	Arrival Time	CT Pneumonia Triage
Vloet, Sanne	98726534	CT	CT Chest W/O Contrast	20 min ago	Suspected Pneumonia
Li, Angel	26378940	CT	CT Chest W/O Contrast	9 min ago	Suspected Pneumonia
Glazer, Aaron	87634538	X-Ray	CXR ABD	3 hr ago	
Miller, Sophia	28394751	X-Ray	CXR Chest	30 min ago	
Rice, Jessica	12361524	CT	CT Head W/O Contrast	15 min ago	

Worklist incorporating CT Pneumonia Triage results

***This product is still in development and has not been reviewed by the FDA. This product is not for sale in the US.**