

**510(k) SUMMARY**  
**KEYA MEDICAL'S DEEPVESSEL PLAQUE**  
**K260497**

**Submitter**

Keya Medical Technology Co., Ltd  
No.1502-1,1502-2,1502-3, Floor 15, Building 1 , Yard No.1, Ronghua South Road  
Beijing Economic-Technological Development Area, Beijing, P.R. China  
Contact Person: Ning Li  
Phone: 011-86-10-87531587  
E-mail: ningli@keyamedical.com

**Date of Preparation:** June 4, 2026

**Name of Device:** DEEPVESSEL Plaque

**Classification Name:** Radiological Image Processing System

**Regulatory Class:** Class II

**Product Code:** QIH

**Regulation:** 21 CFR §892.2050

**Predicate Device:** Autoplaque 3.0 (K212758, Cedars-Sinai Medical Center)

**Device Description**

DEEPVESSEL Plaque ("DVPlaque", the subject device) is a web-based image analysis application for the purpose of characterizing and quantifying plaque and stenosis in coronary arteries based on previously acquired CCTA DICOM data.

DVPlaque applies deep learning neural networks to segment coronary arteries, lumen and vessel walls, which will be reviewed and edited by trained and qualified analysts, if needed. Plaque, stenosis, and vessel measurements are outputted based on the analyst-editable segmentation, utilizing simple rule-based mathematical calculation.

DVPlaque provides a visualization of DVPlaque analysis in two reports — Cardisight Overview Report and Plaque Analysis Report. In Cardisight Overview Report, location and CAD-RADS stenosis category of the plaque are visually over-viewed in a simplified 2D model if stenosis and plaque exist.

In Plaque Analysis Report, quantitative characteristics of the plaque, including total plaque volume, calcified plaque (CP) volume, non-calcified plaque (NCP) volume, low-density plaque (LDP) volume, and plaque burden (defined as percent plaque volume and calculated as total plaque volume/vessel volume) are outputted both on main vessel branch level and on patient level.

Reports of DVPlaque are provided to clinicians to enable them to assess the presence and extent of coronary plaque, aiding in determining treatment paths. These reports are not intended to be the final report used in patient diagnosis and treatment, and should be reviewed with other clinical information. Clinicians may request a re-analysis of the CCTA DICOM data if they do not agree with the report analyses.

### **Intended Use / Indications for Use**

DEEPVESSEL Plaque is a coronary analysis software developed for the analysis of Computed Tomography Angiography (CTA) DICOM data for patients aged 22 and above with suspected or known coronary artery disease (CAD). It provides presence and extent of coronary plaques and stenosis in patients who underwent CCTA for evaluation of suspected or known CAD.

Users should be aware that certain views make use of interpolated data. This is data that is created by software based on the original data set. Interpolated data may give the appearance of healthy tissue in situations where pathology that is near or smaller than the scanning resolution may be present.

The results of the analysis are provided to support qualified clinicians to aid in the evaluation and assessment of coronary arteries. DEEPVESSEL Plaque results are intended to be used by qualified clinicians in conjunction with the patient's clinical history, symptoms, and other diagnostic tests, as well as the clinician's professional judgment. The software is not intended to replace the skill and judgement of a qualified medical practitioner and should only be used by people who have been appropriately trained to use the software.

### ***Comparison to the Predicate Device***

The overall intended use for DEEPVESSEL Plaque is the same as the predicate device: both devices are intended to process radiological images and provide measurements of anatomical components in those images. Furthermore, the subject device's indications for use statement is similar to the predicate device. Both devices are indicated to provide an optimized non-invasive application to analyze coronary anatomy and pathology from CT angiographic images, to provide a post-processing tool for viewing and analyzing cardiac CT data for determining the presence and extent of coronary plaques and luminal stenosis, and to provide analysis results to aid in the determination of treatment paths by the healthcare professional, in conjunction with other patient data. Additionally, both devices are indicated for use by trained operators.

The differences in indications for use between the subject device and predicate device do not cause the subject device to have a different overall intended use. In particular, the location of use for the software operator differs, as DEEPVESSEL Plaque is indicated as a web-based application that can be operated off-site with imaging data transmitted via web-based software, while the predicate device is indicated as a workstation-based desktop device for use at clinical sites. However, for both devices the final analysis results are indicated to be provided to a healthcare professional for use at a clinical location. Therefore, the indicated recipients and locations of application of the final device output are the same for both devices. Of note, the location of use and workflow of the DEEPVESSEL Plaque software is the same as the reference device – i.e., trained operators located in regional data centers use the software device to perform an analysis and generate a report that is subsequently provided to a healthcare professional end user.

### **Comparison of Technological Characteristics**

The technology of DEEPVESSEL Plaque is similar to the predicate device (see **Table 1**). Both devices analyze CT angiographic images for the presence and extent of coronary plaque and luminal stenosis. In addition, they both require trained operators to perform the assessment. Both devices utilize deep learning-based vessel, plaque and lumen segmentation, which is reviewed and can be edited, if necessary, by the operators. Both devices provide equivalent quantitative measurements in a report to the healthcare professional after completing the analysis. These reports are intended to be reviewed by the healthcare professional and used in conjunction with other patient information, such as the original CT images, clinical history, symptoms, clinical risk factors, results of other diagnostic tests, and the clinical judgement of the ordering prescriber. Both devices have demonstrated agreement between clinical expert readers and software measurements in clinical performance testing and have demonstrated appropriate inter-operator and intra-operator agreement.

There are two key technological differences between the subject device and the predicate. First, the location of the DICOM image data processing is different: The DEEPVESSEL Plaque Cloud software is web-based and receives DICOM images via transmittal by a secure picture archiving and communication system (PACS), while the AutoPlaque predicate software resides on a commercial computer platform at the clinical site and receives DICOM image data via loading from the local computer hard drive. Second, the location of the trained software operator may be different: The subject device operator can be at a remote site, while the predicate device operator is at the clinical site on a desktop computer. For both devices, the software operators are trained to review anatomy, define and correct the vessels of interest and generate a report. Both devices allow the end user to revise the analysis as needed, but the approach differs: with AutoPlaque, the user can directly repeat or modify the analysis on-site; with DEEPVESSEL Plaque, the user can request the remote operator to perform a re-analysis.

Therefore, the minor technological differences between the subject and predicate device do not raise different questions of safety or effectiveness.

**Table 1. Technology Comparison to Predicate**

| <b>Parameter</b>                   | <b>Subject Device:<br/>DEEPVESSEL Plaque</b>                                                                                                           | <b>Proposed Predicate<br/>Device:<br/>Autoplaque 3.0<br/>(K212758)</b>  | <b>Comments</b>                                                                                                                                                                                 |
|------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Computer Operating System          | Client-Server Google Chrome Application                                                                                                                | Windows OS<br>Mac OS                                                    | Substantially similar. The difference in platform does not impact safety or effectiveness of the device.                                                                                        |
| Stand-alone software               | Yes                                                                                                                                                    | Yes                                                                     | Same                                                                                                                                                                                            |
| DICOM Compliance                   | DICOM 3.0 Compliant (or higher)                                                                                                                        | DICOM 3.1 Compliant (or higher)                                         | Substantially similar. The difference in platform does not impact safety or effectiveness of the device.                                                                                        |
| 2D Imaging                         | Review of coronary vessels in 2D Multi-Planar Reconstruction (MPR), Curved planar reformation (CPR), and Straightened Curved Planar Reformation (sCPR) | Review of coronary vessels in 2D MPR, curved MPR, and straightened view | Same. The difference is in the nomenclature only.                                                                                                                                               |
| 2D Measurement                     | 2D measurement tools of vessel diameter and contour.                                                                                                   | 2D measurement tools of vessel diameter and contour                     | Same                                                                                                                                                                                            |
| 3D Imaging                         | Yes                                                                                                                                                    | Review of structures in 3D                                              | Same. The difference is in the nomenclature only.                                                                                                                                               |
| Maximum intensity projection (MIP) | No                                                                                                                                                     | Yes                                                                     | This visualization method is not necessary analyzing the presence and extent of coronary plaque and luminal stenosis. The difference in platform does not impact safety or effectiveness of the |

| Parameter                                                                            | Subject Device:<br>DEEPVESSEL Plaque | Proposed Predicate<br>Device:<br>Autoplaque 3.0<br>(K212758) | Comments                                                                                                |
|--------------------------------------------------------------------------------------|--------------------------------------|--------------------------------------------------------------|---------------------------------------------------------------------------------------------------------|
|                                                                                      |                                      |                                                              | device.                                                                                                 |
| Multiplanar reformatting (MPR): MPR with oblique slicing and variable slab thickness | No                                   | Yes                                                          | This is not necessary analyzing the presence and extent of coronary plaque and luminal stenosis.        |
| Volume rendering based on centerlines                                                | Yes                                  | No                                                           | A preprocessing method that modifies the lumen boundary, with no impact on product efficacy or safety.  |
| SCCT Segment Naming                                                                  | Yes (user edit)                      | Yes (user edit)                                              | Same                                                                                                    |
| Automatic/Semi-automatic lumen boundary determination                                | Yes (user edit)                      | Yes (user edit)                                              | Same                                                                                                    |
| Visualize plaque information                                                         | Yes (user edit)                      | Yes (user edit)                                              | Same                                                                                                    |
| Ideal Model Adjustment                                                               | Yes                                  | No                                                           | A graphical representation based on plaque information that does not affect product efficacy or safety. |
| Graphic and text results                                                             | Yes                                  | Yes                                                          | Same                                                                                                    |
| <b>Quantitative Measurements</b>                                                     |                                      |                                                              |                                                                                                         |
| Vessel, plaque, and lumen segmentation                                               | Yes                                  | Yes                                                          | Same                                                                                                    |
| Total Plaque Volume                                                                  | Yes                                  | Yes                                                          | Same                                                                                                    |
| Calcified Plaque Volume                                                              | Yes                                  | Yes                                                          | Same                                                                                                    |
| NCP volume: noncalcified plaque volume                                               | Yes                                  | Yes                                                          | Same                                                                                                    |

| Parameter                                                                                  | Subject Device:<br>DEEPVESSEL Plaque | Proposed Predicate<br>Device:<br>Autoplaque 3.0<br>(K212758) | Comments                                                               |
|--------------------------------------------------------------------------------------------|--------------------------------------|--------------------------------------------------------------|------------------------------------------------------------------------|
| LD-NCP volume: low-density noncalcified plaque volume                                      | Yes                                  | Yes                                                          | Same                                                                   |
| Vessel Volume                                                                              | No                                   | Yes                                                          | This difference does not impact safety or effectiveness of the device. |
| NCP burden: noncalcified plaque volume / analyzed vessel volume                            | Yes                                  | Yes                                                          | Same                                                                   |
| LD-NCP burden: low-density noncalcified plaque volume / analyzed vessel volume             | Yes                                  | Yes                                                          | Same                                                                   |
| CP burden: calcified plaque volume / analyzed vessel volume                                | Yes                                  | Yes                                                          | Same                                                                   |
| Total plaque burden: total plaque volume / analyzed vessel volume                          | Yes                                  | Yes                                                          | Same                                                                   |
| Plaque composition NCP: noncalcified plaque composition (NCP volume / total plaque volume) | No                                   | Yes                                                          | This difference does not impact safety or effectiveness of the device. |
| Plaque composition CP: calcified plaque composition (CP volume / total plaque volume)      | No                                   | Yes                                                          | This difference does not impact safety or effectiveness of the device. |
| Plaque composition LDNCP: low-density noncalcified plaque                                  | No                                   | Yes                                                          | This difference does not impact safety or effectiveness of the device. |

| Parameter                                                                                                            | Subject Device:<br>DEEPVESSEL Plaque | Proposed Predicate<br>Device:<br>Autoplaque 3.0<br>(K212758) | Comments                                                                     |
|----------------------------------------------------------------------------------------------------------------------|--------------------------------------|--------------------------------------------------------------|------------------------------------------------------------------------------|
| composition (LD-NCP<br>volume / NCP volume)                                                                          |                                      |                                                              |                                                                              |
| Diameter stenosis:<br>Maximum diameter<br>stenosis, with respect to<br>proximal and distal<br>references             | Yes                                  | Yes                                                          | Same                                                                         |
| QCAD: Maximum<br>diameter stenosis                                                                                   | Yes                                  | Yes                                                          | Same                                                                         |
| Remodeling index: ratio of<br>maximum vessel area /<br>proximal and distal<br>references                             | No                                   | Yes                                                          | This difference does not<br>impact safety or<br>effectiveness of the device. |
| Area stenosis: maximum<br>area stenosis, with<br>respect to proximal and<br>distal references                        | No                                   | Yes                                                          | This difference does not<br>impact safety or<br>effectiveness of the device. |
| Plaque length: diseased<br>vessel length                                                                             | No                                   | Yes                                                          | This difference does not<br>impact safety or<br>effectiveness of the device. |
| Contrast density<br>difference: maximum<br>difference in contrast<br>density over lesion with<br>respect to proximal | No                                   | Yes                                                          | This difference does not<br>impact safety or effectiveness<br>of the device. |
| MLD: minimal luminal<br>dimeter over lesion                                                                          | Yes                                  | Yes                                                          | Same                                                                         |
| MLA: minimum luminal<br>area over lesion                                                                             | No                                   | Yes                                                          | This difference does not<br>impact safety or<br>effectiveness of the device. |

| Parameter                                                                                                        | Subject Device:<br>DEEPVESSEL Plaque | Proposed Predicate<br>Device:<br>Autoplaque 3.0<br>(K212758) | Comments                                                                     |
|------------------------------------------------------------------------------------------------------------------|--------------------------------------|--------------------------------------------------------------|------------------------------------------------------------------------------|
| Vessel profile: Area,<br>maximum diameter,<br>minimum diameter<br>measured from selected<br>vessel cross section | No                                   | Yes                                                          | This difference does not<br>impact safety or<br>effectiveness of the device. |
| Lumen profile: Area,<br>maximum diameter,<br>minimum diameter<br>measured from selected<br>lumen cross section   | No                                   | Yes                                                          | This difference does not<br>impact safety or<br>effectiveness of the device. |

## Performance Testing

The 510(k) submission provided performance data to establish the substantial equivalence of DVPlaque to the predicate device. A summary of these performance tests is provided below.

### ***Software Verification and Validation:***

Software verification and validation testing were conducted to demonstrate that the subject device meets specifications and works as intended. This included algorithm unit testing, algorithm integration testing, full system testing, and service portal integration testing.

### ***Reproducibility/Repeatability Evaluations:***

Reproducibility & Repeatability (R&R) testing was performed on a group of CT scans with diverse disease conditions to evaluate the variation of repeated analyses of DEEPVESSEL Plaque with different image analysts (reproducibility) at different days with a washout-period in between to avoid memory effects (repeatability). The R&R dataset represented included varying plaque volumes, plaque compositions, and stenosis grades 0 through 4. The R&R cases included a total of 23 subjects with a mean age of 62.2 years. Cases came from the same three U.S. clinical institutions used for validation, and scans represented GE, Arineta, Toshiba, and Siemens CT manufacturers. Testing results met the pre-specified variability metric threshold and thus demonstrated acceptable performance. See **Table 2** below.

**Table 2. Summary of Repeatability and Reproducibility Testing**

| Endpoint            | Factor  | Acceptance Criterion           | Pass/Fail |
|---------------------|---------|--------------------------------|-----------|
| Total Plaque Volume | Day     | %CV due to Day $\leq$ 5.0%     | Pass      |
| Total Plaque Volume | Analyst | %CV due to Analyst $\leq$ 5.0% | Pass      |
| Vessel Volume       | Day     | %CV due to Day $\leq$ 5.0%     | Pass      |
| Vessel Volume       | Analyst | %CV due to Analyst $\leq$ 5.0% | Pass      |

***Clinical Validation Study:***

A retrospective, multi-center clinical validation study was conducted to evaluate the performance of the DEEPVESSEL Plaque software for the detection, quantification, and characterization of coronary atherosclerotic plaque and stenosis using coronary computed tomography angiography (CCTA).

The clinical validation study included a total of 147 CCTA datasets collected from three independent U.S. clinical sites.

The 147 subjects included in the study ranged from 22 to 88 years old and included 91 men and 56 women. Ethnicity information was collected with 72 White (49.0%), 35 Hispanic (23.8%), 22 Asian (15.0%), 4 Black (2.7%), 1 others (0.7%) and unknown 13 (8.8%). The validation cohort also represented a broad range of coronary disease severity, with patient-level CAD-RADS categories spanning CAD-RADS 0-5, including 61 subjects (41.5%) with clinically significant stenosis (CAD-RADS  $\geq$ 3). In addition, all major plaque subtypes were represented in the dataset, including calcified plaque (130 subjects; 88.4%), non-calcified plaque (134 subjects; 91.2%), and low-density plaque (133 subjects; 90.5%).

Each subject contributed one CCTA image dataset. Accordingly, vessel-level analyses were performed on 441 coronary vessel territories, consisting of the left main plus left anterior descending artery (LM+LAD), left circumflex artery (LCX), and right coronary artery (RCA) for each subject, and plaque-level analyses evaluated 711 DVPlaque-detected plaques against expert-reader ground-truth plaques in the primary plaque-detection analysis. CCTA studies were acquired using  $\geq$ 64-detector row CT scanners from multiple manufacturers, including GE, Toshiba, Arineta, and Siemens, representing a range of scanner models and imaging parameters.

Performance of DEEPVESSEL Plaque was evaluated against a reference standard established by consensus of three expert readers. Two expert readers independently reviewed each case and

provided measurements of plaque volume, plaque composition, vessel and lumen volumes, and CAD-RADS stenosis category. The third senior expert reader adjudicated discrepancies according to predefined study rules to establish the final consensus truth. This consensus served as the ground truth for all analyses.

The pivotal clinical validation dataset was independent from the development dataset. The development and validation datasets were separated by site, with development data comprising CCTA images from 1,200 patients acquired from five medical centers in China and one medical center in Europe, while the pivotal validation dataset was collected from three U.S. clinical sites. The protocol also excluded data previously analyzed by DVPlaque. DVPlaque analysts and expert readers were blinded to each other's results during validation.

Clinical performance evaluation demonstrated strong agreement between DEEPVESSEL Plaque and the reference standard for quantitative plaque measurements and stenosis assessment. Correlation analyses showed high agreement for total plaque volume, calcified and non-calcified plaque volumes, vessel volume, and lumen volume at both the patient and vessel levels. Agreement for CAD-RADS stenosis category demonstrated substantial concordance with expert assessment. The results are summarized in **Table 3**. All performance testing results met pre-defined acceptance criteria.

Subgroup analyses demonstrated consistent performance across coronary vessel territories, stenosis severity categories, scanner manufacturers, imaging parameters, and patient demographics, indicating robust and generalizable performance across a clinically representative population.

**Table 3. Key Endpoints and Performance of DVPlaque**

| Performance Metric<br>(Per-Patient) | Statistic                       | Results | Pass/Fail |
|-------------------------------------|---------------------------------|---------|-----------|
| Total Plaque Volume                 | Pearson Correlation Coefficient | 0.976   | Pass      |
| Calcified Plaque Volume             | Pearson Correlation Coefficient | 0.994   | Pass      |
| Non-calcified Plaque Volume         | Pearson Correlation Coefficient | 0.945   | Pass      |
| Vessel Volume                       | Pearson Correlation Coefficient | 0.985   | Pass      |
| Plaque Burden                       | Pearson Correlation Coefficient | 0.978   | Pass      |
| Stenosis Category                   | Weighted Kappa                  | 0.810   | Pass      |
| Vessel Segmentation                 | Dice coefficient                | 0.888   | Pass      |
| Lumen Segmentation                  | Dice coefficient                | 0.887   | Pass      |

## **Conclusion**

The information provided above supports that DEEPVESSEL Plaque is as safe and effective as the predicate device. The device has the same intended use as the predicate device, and minor differences in indications for use do not alter the intended use of the device and do not affect its safety and effectiveness. In addition, the minor technological differences between the devices do not raise different questions of safety or effectiveness, and they are addressed through software verification and validation testing and a clinical validation study. Therefore, DEEPVESSEL Plaque is substantially equivalent to the predicate device.